

**Cold Filter Plugging Point (CFPP) Standard acc to ASTM D6371 and EN 116 (winter diesel)**

Art. ID LM24-OP-4012527  
Unit 250 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cold filter plugging po int (CFPP)		(-29,3 ± 0,3)	°C	ASTM D6371 / EN 116	