

## **Viscosity Standard acc to ASTM D445 and EN ISO 3104 (S3)**

Art. ID LM24-OP-4012536  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

### Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Kinematic Viscosity at 40 °C		2,917 ± 0,003	mm <sup>2</sup> /s	ASTM D445 / EN ISO31 04	
	Dynamic Viscosity		2,372 ± 0,003	mPa.s	calculatio n	