

Custom Solution - 11 component in HNO3 5%

Art. ID LM24-CUS-71578

Unit 125 mL

Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	mg/L		
	Cadmium (Cd)	[7440-43-9]	2	mg/L		
	Chromium (Cr)	[7440-47-3]	10	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		
	Lead (Pb)	[7439-92-1]	30	mg/L		
	Selenium (Se)	[7782-49-2]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	20	mg/L		