

Eltra Computer System - Dell OptiPlex or similar - 10th Generation Intel® Core™ i5-10500 (6-Core, 12MB Cache, 3.1GHz to 4.5GHz, 65W) - Windows 10 Pro 64bit English - 8GB, 1X8GB, DDR4 non-ECC Memory - M.2 256GB PCIe NVMe Class 35 Solid State Drive - A 23.8" monitor enhances your daily workflow. Featuring a height adjustable stand, Full HD resolution and integrated - speakers in a space-saving design.

Art. ID AR-71015
Unit each