

## **Bilirubin**

Art. ID                      NIST-916b  
Unit                        100 mg  
Deliverydetails            Dry ice shipment /not restricted

### **Description**

This Standard Reference Material (SRM) is intended for use in the calibration and standardization of measurement procedures for the determination of bilirubin in clinical samples, and for assignment of bilirubin quantity values to in-house control materials. This SRM comprises neat unconjugated bilirubin that is certified as a chemical substance of known purity. A unit of SRM 916b consists of 100 mg of bilirubin. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Mass Fraction	Bilirubin	[635-65-4]	94,7 ± 0,7	%		