

Custom Solution - 8 components 100 mg/L each in HNO3 3%, HCl 9%

Art. ID LM24-CUS-71590
Unit 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Tungsten (W)	[7440-33-7]	100	mg/L		
	Zirconium (Zr)	[7440-67-7]	100	mg/L		