

Custom Solution - 9 components 100 mg/L each in HNO3 3%, HCl 9%

Art. ID LM24-CUS-71591
Unit 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Gallium (Ga)	[7440-55-3]	100	mg/L		
	Indium (In)	[7440-74-6]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		