

Custom Solution - 9 components 100 mg/L each in HNO3 3%, HCl 9%

Art. ID LM24-CUS-71591
Unit 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Bismuth (Bi) | [7440-69-9] | 100 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 100 | mg/L | | |
| | Gallium (Ga) | [7440-55-3] | 100 | mg/L | | |
| | Indium (In) | [7440-74-6] | 100 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 100 | mg/L | | |
| | Tin (Sn) | [7440-31-5] | 100 | mg/L | | |
| | Tellurium (Te) | [13494-80-9] | 100 | mg/L | | |
| | Titanium (Ti) | [7440-32-6] | 100 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 100 | mg/L | | |