

Custom Solution - 8 components 100 mg/L each in HNO3 5%

Art. ID LM24-CUS-71592
Unit 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	100	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		