

Tolterodine Dimer
(2,2'-((isopropylazanediy)bis(1-phenylpropane-3,1-diyl))bis(4-methylphenol))

Art. ID USP-1A01620

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|---------------|---------------------|------|--------|--------|
| | Tolterodine Dimer | [854306-72-2] | | | | |