

Custom Solution - ICP-MS Tuning solution 25 components in HNO3 2%

Art. ID LM24-CUS-71607
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	µg/L		
	Silver (Ag)	[7440-22-4]	6	µg/L		
	Barium (Ba)	[7440-39-3]	4	µg/L		
	Beryllium (Be)	[7440-41-7]	35	µg/L		
	Bismuth (Bi)	[7440-69-9]	3	µg/L		
	Cerium (Ce)	[7440-45-1]	3	µg/L		
	Cobalt (Co)	[7440-48-4]	8	µg/L		
	Caesium (Cs)	[7440-46-2]	3	µg/L		
	Copper (Cu)	[7440-50-8]	15	µg/L		
	Gallium (Ga)	[7440-55-3]	10	µg/L		
	Holmium (Ho)	[7440-60-0]	3	µg/L		
	Indium (In)	[7440-74-6]	3	µg/L		
	Lithium (Li)	[7439-93-2]	8	µg/L		
	Magnesium (Mg)	[7439-95-4]	10	µg/L		
	Manganese (Mn)	[7439-96-5]	6	µg/L		
	Nickel (Ni)	[7440-02-0]	15	µg/L		
	Rhodium (Rh)	[7440-16-6]	3	µg/L		
	Scandium (Sc)	[7440-20-2]	8	µg/L		
	Strontium (Sr)	[7440-24-6]	5	µg/L		
	Tantalum (Ta)	[7440-25-7]	3	µg/L		
	Terbium (Tb)	[7440-27-9]	3	µg/L		
	Thallium (Tl)	[7440-28-0]	4	µg/L		
	Uranium (U)	[7440-61-1]	3	µg/L		
	Yttrium (Y)	[7440-65-5]	3	µg/L		
	Zinc (Zn)	[7440-66-6]	20	µg/L		