

Custom Solution - 4 components 10000 µg/mL each of in Diethylether

Art. ID LM24-CUS-71229-4x1ML
Unit 4 x 1 mL ampoule
Deliverydetails De Minimis UN1155 I E3

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Dichloroethane D4	[17060-07-0]	10000	µg/mL		
	4-Bromofluorobenzene	[460-00-4]	10000	µg/mL		
	Dibromofluoromethane	[1868-53-7]	10000	µg/mL		
	Toluene D8	[2037-26-5]	10000	µg/mL		