

Custom Solution - 4 components 10000 µg/mL each of in Diethylether

Art. ID LM24-CUS-71229-4x1ML
Unit 4 x 1 mL ampoule
Deliverydetails De Minimis UN1155 I E3

Description

Conservation: 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------|--------------|---------------------|-------|--------|--------|
| | 1,2-Dichloroethane D4 | [17060-07-0] | 10000 | µg/mL | | |
| | 4-Bromofluorobenzene | [460-00-4] | 10000 | µg/mL | | |
| | Dibromofluoromethane | [1868-53-7] | 10000 | µg/mL | | |
| | Toluene D8 | [2037-26-5] | 10000 | µg/mL | | |