

Custom Solution - Farbstoffe - Mix 6 je 20 µg/ml in Ethanol

Art. ID LM24-CUS-10235
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1170 III E1

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Direct Red 28 | [573-58-0] | | µg/mL | | |
| | Direct Black 38 | [1937-37-7] | | µg/mL | | |
| | Basic Blue 26 | [2580-56-5] | | µg/mL | | |
| | Solvent Blue 4 | [6786-83-0] | | µg/mL | | |
| | Basic Red 9 | [569-61-9] | | µg/mL | | |
| | Disperse Blue 1 | [2475-45-8] | | µg/mL | | |