

Custom Solution - 13 components each 100 µg/mL in Methanol

Art. ID LM24-CUS-71611
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fenhexamid	[126833-17-8]	100	µg/mL		
	Fluopicolide	[239110-15-7]	100	µg/mL		
	Imidacloprid	[138261-41-3]	100	µg/mL		
	Iprovalicarb	[140923-17-7]	100	µg/mL		
	Lenacil	[2164-08-1]	100	µg/mL		
	Spiroxamine	[118134-30-8]	100	µg/mL		
	Tebuconazole	[107534-96-3]	100	µg/mL		
	Tebufenozide	[112410-23-8]	100	µg/mL		
	Tetraconazole	[112281-77-3]	100	µg/mL		
	Metolachlor	[51218-45-2]	100	µg/mL		
	Molinate	[2212-67-1]	100	µg/mL		
	Oxadiazon	[19666-30-9]	100	µg/mL		
	Metalaxyl-M	[70630-17-0]	100	µg/mL		