

## **Custom Solution - 12 components each 100 µg/mL in Methanol**

Art. ID LM24-CUS-71612  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dichlorvos	[62-73-7]	100	µg/mL		
	Chlorfenvinphos	[470-90-6]	100	µg/mL		
	Chlorpyrifos	[2921-88-2]	100	µg/mL		
	Azoxystrobin	[131860-33-8]	100	µg/mL		
	Flufenacet	[142459-58-3]	100	µg/mL		
	Atrazine	[1912-24-9]	100	µg/mL		
	Diuron	[330-54-1]	100	µg/mL		
	Isoproturon	[34123-59-6]	100	µg/mL		
	Linuron	[330-55-2]	100	µg/mL		
	Simazine	[122-34-9]	100	µg/mL		
	Terbutylazine	[5915-41-3]	100	µg/mL		
	Terbutryn	[886-50-0]	100	µg/mL		