

Custom Solution - 1 component, Cyanides (CN-) 0,5 mg/L in Potassium Hydroxide 0,1%

Art. ID LM24-CUS-71627
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1814 II E2

Description

Conservation 8 months. (4 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Cyanides (CN-) | | 0,5 | mg/L | | |