

**Custom Solution - 1 component, Cyanides (CN-) 0,6 mg/L in Potassium Hydroxide 0,1%**

Art. ID LM24-CUS-71628  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1814 II E2

Description

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cyanides (CN-)		0,6	mg/L		