

**Custom Solution - 16 components 50 mg/L each of in Acetone + Extended Shelf-Life (Option D): + 12 months**

Art. ID LM24-CUS-71635  
Unit 1 x 1,2 mL ampoule  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation 12 months + Extended Shelf-Life (Option D) + 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alachlor	[15972-60-8]	50	mg/L		
	Atrazine	[1912-24-9]	50	mg/L		
	trans-Chlordane	[5103-74-2]	50	mg/L		
	4,4'-DDD	[72-54-8]	50	mg/L		
	4,4'-DDE	[72-55-9]	50	mg/L		
	4,4'-DDT	[50-29-3]	50	mg/L		
	Endrin	[72-20-8]	50	mg/L		
	beta-HCH	[319-85-7]	50	mg/L		
	Aldrin	[309-00-2]	50	mg/L		
	cis-Chlordane	[5103-71-9]	50	mg/L		
	2,4'-DDD	[53-19-0]	50	mg/L		
	2,4'-DDE	[3424-82-6]	50	mg/L		
	2,4'-DDT	[789-02-6]	50	mg/L		
	Dieldrin	[60-57-1]	50	mg/L		
	alpha-HCH	[319-84-6]	50	mg/L		
	gamma-HCH (Lindan)	[58-89-9]	50	mg/L		