

**Custom Solution - 1 component Monolinuron [CAS:1746-81-2] 100 mg/L in Acetone + Extended Shelf-Life (Option D): + 12 months**

Art. ID LM24-CUS-71642  
Unit 5 x 1,2 mL ampoule  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

Conservation 12 months + Extended Shelf-Life (Option D) + 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Monolinuron	[1746-81-2]	100	mg/L		