

Custom Solution - 13 components 10 mg/L each of in HNO3 2%

Art. ID LM24-CUS-71648
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Germanium (Ge)	[7440-56-4]	10	mg/L		
	Hafnium (Hf)	[7440-58-6]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Niobium (Nb)	[7440-03-1]	10	mg/L		
	Silicon (Si)	[7440-21-3]	10	mg/L		
	Silver (Ag)	[7440-22-4]	10	mg/L		
	Tantalum (Ta)	[7440-25-7]	10	mg/L		
	Tellurium (Te)	[13494-80-9]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Titanium (Ti)	[7440-32-6]	10	mg/L		
	Tungsten (W)	[7440-33-7]	10	mg/L		
	Zirconium (Zr)	[7440-67-7]	10	mg/L		