

### **Custom Solution - 8 components each at 1000 mg/L in Methanol**

Art. ID LM24-CUS-71657  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aniline	[62-53-3]	1000	mg/L		
	o-Toluidine	[95-53-4]	1000	mg/L		
	p-Toluidine	[106-49-0]	1000	mg/L		
	Nitrobenzene	[98-95-3]	1000	mg/L		
	2-Nitrotoluene	[88-72-2]	1000	mg/L		
	3-Nitrotoluene	[99-08-1]	1000	mg/L		
	4-Nitrotoluene	[99-99-0]	1000	mg/L		
	m-Toluidine	[108-44-1]	1000	mg/L		