

Custom Solution - N(NH₄⁺) 2000 mg/L in Water

Art. ID LM24-CUS-71670
Unit 1 x 100 mL bottle
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 24 months. (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|----------------------------------|--|------------|---------------------|------|--------|--------|
| N(NH ₄ ⁺) | N-(Ammonium (NH ₄ ⁺)) | | 2000 | mg/L | | |