

Iron Base, Steel - Stainless, CRM, 304 H, , Disc ~47 mm x ~30 mm

Art. ID VS-LG61
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,033 | % | | |
| | Carbon (C) | [7440-44-0] | 0,307 | % | | |
| | Chromium (Cr) | [7440-47-3] | 18,8 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,065 | % | | |
| | Manganese (Mn) | [7439-96-5] | 1,51 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 1,05 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,47 | % | | |
| | Nickel (Ni) | [7440-02-0] | 9,18 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0133 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0107 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,83 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,4 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,197 | % | | |
| | Tungsten (W) | [7440-33-7] | 1,11 | % | | |