

PFAS EF-28 Native Mixture in Methanol

Art. ID CIL-ES-5636-A
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1H,1H,2H,2H-Perfluorodecanesulfonate (8:2 FTS), sodium salt	[27619-96-1]	1000	ng/mL		
	1H,1H,2H,2H-Perfluorodecanesulfonate (10:2 FTS), sodium salt	[108026-35-3]	1000	ng/mL		
	1H,1H,2H,2H-Perfluorooxanesulfonate (4:2 FTS), sodium salt	[27619-93-8]	1000	ng/mL		
	2-(Perfluorohexyl)ethanesulfonic Acid Sodium Salt	[27619-94-9]	1000	ng/mL		
	N-Ethyl-N-(perfluoro-1-octanesulfonyl) Glycine	[2991-50-6]	1000	ng/mL		
	N-Methylperfluoro-1-octanesulfonamidoacetic Acid	[2355-31-9]	1000	ng/mL		
	Bis[2-(perfluorooctyl)ethyl] Phosphate	[678-41-1]	1000	ng/mL		
	N-Methylperfluorooctanesulfonamide	[31506-32-8]	1000	ng/mL		
	Perfluoro-n-butane sulfonic acid potassium salt	[29420-49-3]	1000	ng/mL		
	Perfluorobutyric Acid, Sodium Salt	[2218-54-4]	1000	ng/mL		
	PERFLUORODECANESULFONATE, POTASSIUM SALT	[2806-16-8]	1000	ng/mL		
	Sodium Perfluorodecanoate	[3830-45-3]	1000	ng/mL		

te			
Perfluorododecanoic Acid, Sodium Salt	[307-67-5]	1000	ng/mL
Sodium perfluoro-1-heptanesulfonate	[21934-50-9]	1000	ng/mL
Perfluoro-n-heptanoic acid	[375-85-9]	1000	ng/mL
Perfluorohexadecanoic acid	[67905-19-5]	1000	ng/mL
Perfluoro-n-hexane sulfonic acid potassium salt	[3871-99-6]	1000	ng/mL
Perfluorohexanoic acid (PFHxA), sodium salt	[2923-26-4]	1000	ng/mL
Perfluoro-n-nonanoic acid	[375-95-1]	1000	ng/mL
Perfluorooctadecanoic acid	[16517-11-6]	1000	ng/mL
Perfluorooctane sulfonamide PFOSA	[754-91-6]	1000	ng/mL
Perfluorooctanoic acid	[335-67-1]	1000	ng/mL
Perfluoropentanesulfonate (PFPeS), sodium salt	[630402-22-1]	1000	ng/mL
Perfluoropentanoic acid (PFPeA), sodium salt	[2706-89-0]	1000	ng/mL
Perfluoro-n-tetradecanoic acid	[376-06-7]	1000	ng/mL
Perfluoro-n-tridecanoic acid	[72629-94-8]	1000	ng/mL
Perfluoroundecanoic acid (PFUnA), sodium salt	[60871-96-7]	1000	ng/mL
Perfluorooctane sulfonic acid sodium salt (PFO S)	[4021-47-0]	1000	ng/mL