

Custom solution - 20 components in 2% HNO₃

Art. ID LM24-CUS-71687
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months. (4 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1	mg/kg		
	Barium (Ba)	[7440-39-3]	1	mg/kg		
	Chromium (Cr)	[7440-47-3]	3,6	mg/kg		
	Copper (Cu)	[7440-50-8]	5	mg/kg		
	Lead (Pb)	[7439-92-1]	1,5	mg/kg		
	Manganese (Mn)	[7439-96-5]	0,6	mg/kg		
	Titanium (Ti)	[7440-32-6]	0,02	mg/kg		
	Antimony (Sb)	[7440-36-0]	0,04	mg/kg		
	Cadmium (Cd)	[7440-43-9]	0,1	mg/kg		
	Cobalt (Co)	[7440-48-4]	0,05	mg/kg		
	Iron (Fe)	[7439-89-6]	48,01	mg/kg		
	Lithium (Li)	[7439-93-2]	0,6	mg/kg		
	Nickel (Ni)	[7440-02-0]	0,02	mg/kg		
	Zinc (Zn)	[7440-66-6]	5	mg/kg		
	Arsenic (As)	[7440-38-2]	0,02	mg/kg		
	Lanthanum (La)	[7439-91-0]	0,0125	mg/kg		
	Europium (Eu)	[7440-53-1]	0,0125	mg/kg		
	Gadolinium (Gd)	[7440-54-2]	0,0125	mg/kg		
	Terbium (Tb)	[7440-27-9]	0,0125	mg/kg		
	Molybdenum (Mo)	[7439-98-7]	0,12	mg/kg		