

Custom solution - Siloxane - Mix 2, each at 100 µg/mL in Acetone

Art. ID LM24-CUS-10237
Unit 1 x 1 mL
Deliverydetails De Minimis UN1090 II E2

Description

Haltbarkeit 18 Monate

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Octamethyltrisiloxane	[107-51-7]	100	µg/mL		
	Octamethylcyclotetrasiloxane	[556-67-2]	100	µg/mL		