

Custom Solution - 12 components in HNO₃ 10%

Art. ID LM24-CUS-71700
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1000	mg/L		
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Barium (Ba)	[7440-39-3]	1000	mg/L		
	Bismuth (Bi)	[7440-69-9]	1000	mg/L		
	Cobalt (Co)	[7440-48-4]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Manganese (Mn)	[7439-96-5]	1000	mg/L		
	Arsenic (As)	[7440-38-2]	500	mg/L		
	Beryllium (Be)	[7440-41-7]	500	mg/L		
	Chromium (Cr)	[7440-47-3]	500	mg/L		
	Copper (Cu)	[7440-50-8]	500	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		