

Custom Solution - 13 components in HNO3 10%

Art. ID LM24-CUS-71701
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Selenium (Se)	[7782-49-2]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Titanium (Ti)	[7440-32-6]	1000	mg/L		
	Thallium (Tl)	[7440-28-0]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Lithium (Li)	[7439-93-2]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Lead (Pb)	[7439-92-1]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	500	mg/L		
	Tellurium (Te)	[13494-80-9]	500	mg/L		