

Custom solution - 8 components 100 mg/L each of in HNO3 5%

Art. ID LM24-CUS-71711
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Strontium (Sr)	[7440-24-6]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		