

Custom Solution - 3 component each 40 mg/mL in Water

Art. ID LM24-CUS-71723
Unit 1 x 100 mL bottle
Deliverydetails No Dangerous Good /not restricted

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Thiocyanate (SCN-)		40	mg/L		
	Thiosulphate (S ₂ O ₃ ²⁻)		40	mg/L		
	Iodide (I-)		40	mg/L		