

Poly(ethyleneimine) solution, average Mn ~60,000 by GPC, average Mw ~750,000 by LS, 50 wt. % in H2O (ALDRICH)

Art. ID SAF-181978-18KG
Unit 18 kg
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3082 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Poly(ethyleneimine)	[9002-98-6]	50	wt %		