

41L40 MOD Low Alloy Steel Disc 41 mm D x 40mm

Art. ID BS-SU 41L40
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,02 | % | | |
| | Arsenic (As) | [7440-38-2] | <0,05 | % | | |
| | Boron (B) | [7440-42-8] | <0,005 | % | | |
| | Carbon (C) | [7440-44-0] | 0,39 | % | | |
| | Calcium (Ca) | [7440-70-2] | <0,005 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,008 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,97 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,13 | % | | |
| | Iron (Fe) | [7439-89-6] | 96,6 | % | | |
| | Magnesium (Mg) | [7439-95-4] | <0,005 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,9 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,2 | % | | |
| | Nitrogen (N) | [7727-37-9] | <0,05 | % | | |
| | Niobium (Nb) | [7440-03-1] | <0,05 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,25 | % | | |
| | Oxygen (O) | [7782-44-7] | <0,05 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,01 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,14 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,02 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,26 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,01 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,002 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,002 | % | | |
| | Tungsten (W) | [7440-33-7] | <0,05 | % | | |
| | Zirconium (Zr) | [7440-67-7] | <0,05 | % | | |