

Apatite Nano-Pellet, pressed pellet, diameter 10 mm (Standard for solid-state microanalysis)

| | |
|-----------------|-----------------------------------|
| Art. ID | MY-Apatite-NP-LA-ICP-MS-10MM |
| Unit | each (pressed pellet) |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Pellet for LA-ICP-MS application /// The principle behind LA-ICP-MS (Laser Ablation - Inductively Coupled Plasma - Mass Spectrometry) involves a laser beam removing (ablating) material from a sample and analysing its chemical composition in a mass spectrometer