

**Polymer Mount 1 Inch Round, each. (Holder for 10 mm Pellets)**

Art. ID	MY-PM-1IR-LA-ICP-MS
Unit	each (pressed pellet)
Deliverydetails	No Dangerous Good /not restricted

Description

Pellet for LA-ICP-MS application /// The principle behind LA-ICP-MS (Laser Ablation - Inductively Coupled Plasma - Mass Spectrometry) involves a laser beam removing (ablating) material from a sample and analysing its chemical composition in a mass spectrometer