

### **Custom Solution - 13 components 100 mg/L each of in Methanol**

Art. ID LM24-CUS-71750  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Chloroaniline	[95-51-2]	100	mg/L		
	3-Chloroaniline	[108-42-9]	100	mg/L		
	4-Chloroaniline	[106-47-8]	100	mg/L		
	2,4-Dichloroaniline	[554-00-7]	100	mg/L		
	2,5-Dichloroaniline	[95-82-9]	100	mg/L		
	2,3-Dichloroaniline	[608-27-5]	100	mg/L		
	3,4-Dichloroaniline	[95-76-1]	100	mg/L		
	3,5-Dichloroaniline	[626-43-7]	100	mg/L		
	2-Ethoxyaniline	[94-70-2]	100	mg/L		
	2-Nitroaniline	[88-74-4]	100	mg/L		
	3-Nitroaniline	[99-09-2]	100	mg/L		
	N,N-Dimethylaniline	[121-69-7]	100	mg/L		
	2-Aminonaphthalene	[91-59-8]	100	mg/L		