

L-Phenylalanine (3',5'-13C2, 99%, 2,3,3,2',4',6'-D6,98%, 15N, 98%)

Art. ID CIL-CDNLM-12287-PK
Unit on request

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Phenylalanine (3',5'- 13C2, 99%, 2,3,3,2',4', 6'-D6,98%, 15N, 98%)	[63-91-2]/unlab elled compound				