

Custom Solution - 10 components in HNO3 5% (Extended Shelf-Life (Option D) + 12 months)

Art. ID LM24-CUS-71765
Unit 2 x 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Extended Shelf-Life (Option D) + 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Silver (Ag)	[7440-22-4]	50	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Lead (Pb)	[7439-92-1]	500	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Arsenic (As)	[7440-38-2]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	500	mg/L		
	Selenium (Se)	[7782-49-2]	500	mg/L		
	Boron (B)	[7440-42-8]	500	mg/L		
	Thallium (Tl)	[7440-28-0]	500	mg/L		