

Custom Solution - 10 components in HNO3 5% (Extended Shelf-Life (Option D) + 12 months)

Art. ID LM24-CUS-71766
Unit 2 x 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Extended Shelf-Life (Option D) + 12 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Barium (Ba) | [7440-39-3] | 50 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 50 | mg/L | | |
| | Beryllium (Be) | [7440-41-7] | 50 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 50 | mg/L | | |
| | Aluminium (Al) | [7429-90-5] | 500 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 500 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 500 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 500 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 500 | mg/L | | |
| | Vanadium (V) | [7440-62-2] | 500 | mg/L | | |