

Custom Solution - 10 components in HNO3 5% (Extended Shelf-Life (Option D) + 12 months)

Art. ID LM24-CUS-71766
Unit 2 x 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Extended Shelf-Life (Option D) + 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	50	mg/L		
	Strontium (Sr)	[7440-24-6]	50	mg/L		
	Beryllium (Be)	[7440-41-7]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	50	mg/L		
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Cobalt (Co)	[7440-48-4]	500	mg/L		
	Molybdenum (Mo)	[7439-98-7]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Chromium (Cr)	[7440-47-3]	500	mg/L		
	Vanadium (V)	[7440-62-2]	500	mg/L		