

Isotachysterol3 (COLD SHIPMENT REQUIRED)

(S)-4-methyl-3-((E)-2-((1R,7aR)-7a-methyl-1-((R)-6-methylheptan-2-yl)-2,3,5,6,7,7a-hexahydro-1H-inden-4-yl)vinyl)cyclohex-3-en-1-ol)

Art. ID USP-1A02050
Unit 25 mg
Deliverydetails Dry ice shipment UN2811 III E1

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Isotachysterol3	[22350-43-2]				