

CONDUCTIVITY SOLUTION ISO 17034 - (RM) 0,XX μ S/cm 100 ml bottle

Art. ID	CF-CHFSITCON0.XX100RM
Unit	100 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Customized Conductivity standard solutions produced in compliance with UNI CEI EN ISO 17034:2017 /// VALUES BETWEEN 0,61 AND 0,89 μ S/cm (RM) /// Single bottle, Product shelf life 3 months from purchase /// The values given refer to the reference temperature of 25 °C