

CONDUCTIVITY SOLUTION CRM ISO 17034 - XX,00 μ S/cm 100 ml bottle

Art. ID	CF-CHFSITCON000XX100
Unit	100 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Customized Conductivity standard solutions produced in compliance with UNI CEI EN ISO 17034:2017 /// VALUES BETWEEN 5,01 A 29,99 μ S/cm (EXCEPT: 10,00-20,00-25,00 μ S/cm) /// Single bottle, Product shelf life 3 months from purchase /// The values given refer to the reference temperature of 25 °C