

CONDUCTIVITY SOLUTION CRM ISO 17034 - XXXXX,00 μ S/cm PACK 300ml in 6x50 ml bottles

Art. ID	CF-CHFSITCON0XXXXX300
Unit	6 x 50 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Customized Conductivity standard solutions produced in compliance with UNI CEI EN ISO 17034:2017 /// VALUES BETWEEN 12.881 AND 100.000 μ S/cm (EXCEPT: 20.000, 50.000, 75.000 AND 100000 μ S/cm) /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 25 °C