

**CONDUCTIVITY SOLUTION CRM ISO 17034 - XXXXX,00  $\mu$ S/cm PACK 500ml in 2x250 ml bottles**

Art. ID	CF-CHFSITCON0XXXXX500
Unit	2 x 250 mL
Deliverydetails	No Dangerous Good /not restricted

**Description**

Customized Conductivity standard solutions produced in compliance with UNI CEI EN ISO 17034:2017 /// VALUES BETWEEN 12.881 AND 100.000  $\mu$ S/cm (EXCEPT: 20.000, 50.000, 75.000 AND 100000  $\mu$ S/cm) /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 25 °C