

REFRACTOMETRY SOLUTION CRM ISO 17034 - °BX 40,00 / nD 1.39986 PACK 6x4 ml vials

| | |
|-----------------|-----------------------------------|
| Art. ID | CF-CHFSITBX4006F |
| Unit | 6 x 4 mL VIALS |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 20 °C ///
Refractometry solutions (°BX / nD) produced in compliance with UNI CEI EN ISO 17034:2017