

REFRACTOMETRY SOLUTIONS CRM ISO 17034 - °BX 0-5-12-25-40 PACK 15 vials 3 for each value

Art. ID CF-CHFSITBX0405M15F
Unit 5 x (3 x 4 mL VIALS)
Deliverydetails No Dangerous Good /not restricted

Description

Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 20 °C ///
Refractometry solutions (°BX / nD) produced in compliance with UNI CEI EN ISO 17034:2017