

REFRACTOMETRY SOLUTION CRM ISO 17034 - °BX 10,00 PACK 500ml in 2x250 ml bottles

| | |
|-----------------|-----------------------------------|
| Art. ID | CF-CHFSITBX10500 |
| Unit | 2 x 250 mL |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Product shelf life 6 months from purchase /// The values given refer to the reference temperature of 20 °C ///
Refractometry solutions (°BX / nD) produced in compliance with UNI CEI EN ISO 17034:2017