

**REFRACTOMETRY SOLUTION CRM ISO 17034 - °BX 30,00 PACK 500ml in 2x250 ml bottles**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | CF-CHFSITBX30500                  |
| Unit            | 2 x 250 mL                        |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Product shelf life 6 months from purchase /// The values given refer to the reference temperature of 20 °C ///  
Refractometry solutions (°BX / nD) produced in compliance with UNI CEI EN ISO 17034:2017