

**ISO pH ELECTRODES CLEANER PH OILCLEANER PACK 300ml in 6 x50 ml bottles 8  
General electrodes cleaning, fatty residues removal)**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | CF-CHFISOCLOIL300                 |
| Unit            | 6 x 50 mL                         |
| Deliverydetails | No Dangerous Good /not restricted |

Description

pHmeter auxiliary solutions with ISO traceable certificate /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 25 °C